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Tel. +49-(0)2131-1335605

EN User Manual Digital Scale LS

Thank you for purchasing this device. Please read the instructions carefully before using them for the first time and then follow all the instructions.

- Avoid exceeding the maximum weighing range. Overloading will damage the load cell.
- Allow at least 30 seconds for the scale to warm up after turning it on. The internal components must first stabilize, especially before calibration. Without stable operating temperature, measurement errors will occur.
- Keep the device in a clean environment. Dust, dirt, moisture, vibration, drafts, and proximity to other electronic devices can affect reliability and accuracy.
- Handle the device with care. Carefully place items to be weighed on the weighing surface. Continued rough handling will damage the internal sensors.
- Only use the device on a stable, horizontal and vibration-free surface.

Technical data

Modell	Min- / Max capacity	Readout	Adjustment weight
LS 200	0,2g / 200g	0,01g	2 * 100g
LS 500	0,2g / 500g	0,01g	2 * 200g
LS 1000	1g / 1000g	0,05g	2 * 500g
LS 2000H	1g / 2000g	0,05g (< 1kg) 0,1g (>1kg)	2 * 1000g
LS 3000	2g / 3000g	0,1g	1000g + 2000g

Weighing units: g, oz, dwt, ozt, ct, gn

Power supply: 2 x AAA batteries

Operating temperature: 10° to 30°C, up to 80% R.H.

Shutdown: 2 minutes without new measurement

Tolerance: Within $\pm 0.01\%$ of max capacity

I. Operation

- Open the battery compartment on the back of the scale and insert the AAA batteries, matching the polarity.
- Place the device on a stable and level surface. Remove the two pans above the weighing surface (Or flip them up). Press the [ON/OFF] button. Wait until "0" appears on the display.
- Place the object to be weighed in the center of the weighing platform. The weight determined appears in the display. If the weight exceeds the

maximum capacity, the display shows "O_LD". Such an overload must be avoided at all costs and can permanently damage the scale!

- Use the [MODE] button to switch between the different weight units.

II. TARE Weighing

- Turn on scale as described above.
- Place the “tare item” on the platform.
- Press the [TARE] key and wait until “0.00” is displayed.
- Add the “net-weight-item”.

III. Weighing units

- Press the [MODE] repeatedly to change the weight unit until the required weight unit appears

IV. Calibration

NOTE: The scale has been calibrated by the factory. For the calibration, a preset adjustment weight is required. Without such a adjustment weight the calibration is not possible and leads to erroneous values.

Please ask the merchant for the suitable test weights.

- Turn on the scale by pressing the [ON/OFF] without a weight on the scale. Wait until the display shows "0.00".
- Press and hold the [MODE] key for 3 seconds until the display shows “CAL”. Again press [MODE] key.
- After a few seconds the Display flashes "XXX g" on the display (required adjustment weight).
- Place the requested weight on the surface and wait until the scale indicates the value. Depending on the model the scale now will show “PASS” or a new calibration-value. If a new weight is required remove the old weight and place the new weight on the surface, wait until it shows “PASS”
- As soon as you see “PASS” the calibration is completed successfully.

V. Guarantee / Warranty

G&G grants a 2-year guarantee on the function of the scale after purchase.

Excluded from the guarantee are:

- Errors and damage caused by external influences, e.g. caused by:
- Overloading, dropping objects on the scale or the scale itself, extreme heat, vibration, dust, moisture, or contact with corrosive materials (unless the particular model is specifically designed for it).
- Errors and damage caused by improper use and non-observance of the relevant instructions in the operating instructions.

- Errors in scales that have been technically modified without the consent of G&G.

The guarantee begins on the date of purchase (in a shop) or, in the case of ordered goods, on the date of delivery and is valid upon presentation of a document stating the date of purchase. The guarantee is given in such a way that scales that show a fault due to material or production-related defects are repaired or exchanged free of charge at our discretion. A refund of the purchase price is not possible in the event of a warranty claim.

The costs of the return will be borne by G&G in the case of justified complaints (limited to the country to which the delivery was made).

In the event of a defect, contact us by email to receive a free return shipment label.

Your legal statutory rights, in particular warranty rights, are not limited by this guarantee and continue to apply to you.

Contact information:

G&G GmbH
Novesiastr. 31
41564 Kaarst
Germany

Email: info@gandg.de,
Phone +49-(0)2131-1335605